

SUPER SHAKER BED VIBRATOR MODEL SS12V

IMPORTANT– Please read these instructions carefully before use and retain for future reference. See operating instructions.

-- Use only an approved Sonic Alert extension cord model SBE115, or connect directly to unit.

Warning: To prevent fire or shock hazard do not expose this appliance to rain or moisture.

-- No naked flame sources, such as lighted candles should be placed on the apparatus.

-- It shall not be exposed to dripping or splashing and that no objects filled with liquids, such as vases, shall be placed on the apparatus.

-- Do not expose the unit to any sources of water or moisture, do not place any objects with water or filled with liquids of any type near the unit.

OPERATION

When activated by your **Sonic Boom™ Clock** or any clock with matching DC outlet.

Sonic Alert Super Shaker bed vibrator will wake even the heaviest sleeper.

The **Super Shaker Model SS12V** bed vibrator has extra strong vibrations designed especially for a long life. The **Super Shaker** bed vibrator also has a built-in automatic shut-off sensor so it will not overheat even if accidentally left on.

INSTALLATION

The **Super Shaker** bed vibrator is very safe and easy to use.

1. Place the **Super Shaker** bed vibrator between your mattress and box springs or under your pillow.
2. Plug the Super Shaker bed vibrator into the **Sonic Boom™ Clock** or into any clock with a matching DC outlet.
3. The **Super Shaker** bed vibrator will vibrate until you turn off the alarm or hit the snooze button. The Super Shaker bed vibrator will not overheat if you forget to turn off the alarm. If the Super Shaker bed vibrator is left on for more than approximately 10 minutes, the built-in temperature sensor will turn off the vibrator so it does not become too hot. When the Super Shaker bed vibrator cools down, the temperature sensor will automatically reset the vibrator for normal operation.

SONIC ALERT'S LIMITED WARRANTY

The SUPER SHAKER bed vibrator, Model SS12V is warranted against manufacturing defects in materials and workmanship for five (5) years from the date of purchase. Within this period Sonic Alert will repair or replace the SUPER SHAKER bed vibrator without charge for parts or labor. Simply send the SUPER SHAKER bed vibrator (postpaid) and your sales slip as proof of purchase to:

Customerservice@sonicalert.com • 248-577-5400
850 Stephenson Hwy., Suite 323, Troy MI 48083

-



This icon indicates separate disposal for electrical and electronic equipment.

The product is covered by the legal guarantee of conformity as provided by applicable law.